

New Certified Reference Materials for Chromatography



Use of Certified Reference Materials (CRMs)

Quality Assurance

- Ensuring Accuracy of the Mean and Precision
- Ensuring an Analytical Method provides correct Results

Traceability of results to acknowledged standards

Calibration


Validation of Analytical Methods

Interlaboratory Trials

ISO 17025 accredited labs have to use CRMs on regular base



Accuracy, Accuracy of the Mean & Precision






not correct

precise

but

Accuracy = accuracy of the mean + precision



Accuracy of the mean: systematic errors

Precision: accidental errors (by random)

Accuracy: overall deviation



Definition "Certified Reference Material"

Reference material, accompanied by a **certificate**, one or more of whose **property values are certified**, by a **procedure** which **establishes its traceability** to an accurate realization of the unit in which the property values are expressed, and for which each certified value is **accompanied by an uncertainty** to a **stated level of confidence**.

VIM:1995, ISO guide 30/31



National Metrological Institutes (NMIs)

Providing acknowledged standards, examples for NMIs are:

USA: National Institute of Standards and Technology (NIST)



Germany: Federal Institute for Materials Research and Testing (BAM)





UK National Physics Laboratory






Uncertainty







Typical uncertainties in chemical analysis:

Simple sample, metrological level:	0.05-0.4%		
Simple sample (no Matrix), normaler level	0.5-1.5%		
Analysis in Matrix	2-50%		

Problem: Users prefer low uncertainties

Consequence is: many CRM manufacturers calculate (or estimate) the uncertainty too optimistic (ignoring relevant shares).



Traceability

ISO 17025 Section 5.6: Measurement traceability 5.6.3 Reference standards and reference materials

Definition: property of a value to be related by an unbroken chain of comparison measurements to a reasonable normative.

(VIM 6:10 – ref: BIPM, IEC, IFCC, ISO, IUPAC, IUPAP, OIML, International Vocabulary of Basic and Metrological Terms in Metrological

The uncertainty of all uncertainties of all comparison meas 2 2 ty of the reference.

The length of this "traceability chain" is not limited, but: the uncertainty is increasing!

X % Cu ±Y %

ISO 17025 + ISO Guide 34 Accreditation

ISO 17025: "General Requirements for the Competence of Calibration and Testing Laboratories"

- Many labs do have 17025 Accreditation
- Important: Tracebility, uncertainty, education, infrastructure
- Accreditation is bound to a lab and scope

ISO Guide 34: "General Requirements for the Competence of Reference Material Producers"

Produced in double accredited laboratory fulfilling
ISO/IEC 17025 and
ISO Guide 34

- Only few labs do have ISO Guide 34 Accreditation
- Important: Tracebility, uncertainty, homogenity, stability, shelf life
- Accreditation is bound to a lab and scope

ISO 34 + 17025: "Gold Standard" Accreditation for Reference Material manufacturers

- Only very few labs do have both
- Acknowledgement of competence and reliability

ISO/IEC 17025



STS 490

ISO Guide 34 und ISO/IEC 17025



SRMS 001

Organic CRMs: Challenges

- Numberless organic compounds
- Only few international standards available
- •Traceability?




Alternative approach: qNMR

Quantitative NMR (qNMR) is a Relative Primary Method

NMR spectrum of an equimolar mixture of CH₃Cl and CH₃Br

Signal intensity is:

- proportional to the number of protons
- independent of chemical structure (chemical shift is impacted)





Advantages of Quantitative NMR


Traceability independent of chemical identity Low Uncertainties

Expanded uncertainties between 0.1% and 0.5%
 Non-destructive
 Structural information on target compound
 Information on impurities


Jancke H. NMR als primäre Methode. Nachrichten aus Chemie, Technik und Laboratorium, 1998, 46, 722 Jancke, H. NMR Spectroscopy as a Primary Analytical Method, Document 98/02 to the 4th Session of the CCQM, Sèvres 1998



Quantitative NMR (qNMR): Procedure









Content Assignment by qNMR

Internal standard and the reference standard are weighed together into one NMR tube (with solvent)


Intensities of appropriate signals of the reference standard (I_x) and the internal standard (I_{Std}) are used for calculation Content P_x can be calculated by the following equation:

$$P_{x} = \frac{I_{x}}{I_{Std}} \cdot \frac{N_{Std}}{N_{x}} \cdot \frac{M_{x}}{M_{Std}} \cdot \frac{m_{Std}}{m_{x}} P_{Std}$$

- Number of nuclei contributing to the resonance (N_x, N_{Std})
- Molar masses (M_x, M_{Std})
- Initial weight (m_x, m_{Std}) and purity of internal standard (P_{Std})




Tracebility to NIST Standard Reference Materials



Acenaphthene Fluka 05426



Signal allocation and purity: 2D-NMR: H-H COSY



2D-COSY of 3,5-dinitrobenzoic acid. Signals of impurity 3-nitrobenzoic acid detectable. Easy detection at concentration below 0,1%







Traceability [2]:

NIST SRM 350b (Benzoic acid)

Certificate


This certificate is designed in accordance with ISO Guide 31 [1].

Acenaphthene Substance:

Product no .: 05426 BCBD1028V Lot no.: Formula: C₁₂H₁₀ Molecular mass: 154.21

Traceability (2): NIST SRM 350b (Benzolc acid)

Certificate Issue date: October 19, 2010 Oct/2012 Expiry:



Certified value and uncertainty according to ISO Guide 35 [3] and Eurachem/CITAC Guide [4]

Substance		Certified value	Expanded uncertainty
	Acenaphthene	98.5 % (g/g)	0.3 % (g/g)

Minimum sample: 50 mg is recommended as the minimal sample amount. If less material is

used, it is recommended to increase the certified uncertainty by a factor of two

for half of sample and a factor of four for one fourth mg of sample.

Drying instruction: This material does not require drying before use.

Intended use: Use this CRM as calibrant for chromatography or any other analytical

technique.

The CRM should be stored in the original bottle at room-temperature (20-25°C). After use the bottle should be tightly closed and protected from Storage and handling:

excessive moisture and light.

CRM operations:

Certification body:



Page 1 of 4 Sigma-Aldrich Production GmbH, industriestrasse 25, 9471 Buche/ Switzerland SIGMA-ALDRICH

Certified value

98.5 % (g/g)

Expanded uncertainty, $U = k \cdot u_c$ (k = 2)

0.3 % (g/g)



Organic TraceCERT: CRMs for Chromatography

Acids and Esters

Aldehydes and Aldehyde Derivatives

Amines and Amides (e.g. acryl amide)

Amino Acids

Antibiotics / Drugs (including veterinary drugs e.g. antibiotics)

FA / FAME

Natural Products

Pesticides

Phthalates (plasticizers)

Polyaromatic Hydrocarbons (PAHs)

SVHC (semi-volatile hydrocarbons)




Difference Matrix CRM ⇔ Pure Substance CRM

Matrix CRM ⇒ certified impurities

Pure Substance CRM ⇒ certified content



Aflatoxin B1 in peanuts: usage of CRMs



Non-problematic: Weighing, dilution, quant. HPLC

Problematic: Grinding, Extraction, SPE

How to conect the HPLC result with the whole analytical method? Neat compound CRM (Aflatoxin B1) only for the verification of the accuracy of the mean of step 4 (ONLY of step 4)



Usage of Matrix CRMs



Spiking with labelled Standard: assumption, that the labelled substance behaves like a non-labelled and that interaction with matrix is the same.

Matrix CRM: shows absence of interferences, lost of the sample, systematic errors. Shows completeness of extraction & purification.



BCR-264 Defatted Peanut Meal Aflatoxin B1, high level





EUROPEAN COMMISSION

COMMUNITY BUREAU OF REFERENCE - BCR

CERTIFIED REFERENCE MATERIAL **BCR-264**

CERTIFICATE OF ANALYSIS

AFLATOXIN B ₁ IN DEFATTED PEANUT MEAL			
Mass fraction	Aflatoxin B ₁ ¹⁾	Number of accepted sets of data p	
Certified value	Uncertainty 3)	sets of data p	
206 μg/kg	13 μg/kg	11	

- 1) Expressed on material as supplied.
- This value is the unweighted mean of p accepted mean values, independently obtained by p sets of results.
 The uncertainty is taken as the half width of the 95 % confidence interval of the mean defined in (2). It does not include any uncertainty on the value for the molar absorption coefficient of aflatoxin B₁ in chloroform which has been adopted for this certification (see § 7.2 of the report)

When the reference material is used to assess the performance of a method, the user should refer to the recommendations laid down in the last chapter (instructions for use) of the certification report.

DESCRIPTION OF THE MATERIAL

The material is a finely ground defatted peanut meal. It is supplied in units of about 150 g in sachets sealed under vacuum.

The minimum recommended sample intake is 20 g

INSTRUCTIONS FOR USE

Samples should be stored unopened at - 18 °C. To avoid water condensation they should be allowed to warm to room temperature (e.g. overnight) before opening. The contents should be thoroughly mixed before sub-samples are taken. After opening, the material should be used on the same day.

PARTICIPATING LABORATORIES

- Institut für Angewandte Botanik, Hamburg (D)
- Istituto Superiore di Sanità, Roma (I)
- Leatherhead Food R.A., Leatherhead, Surrey (UK)
- LUFA, Kiel (D)
- R.I.V.M., Bilthoven (NL)
- RHM Research and Engineering Ltd., High Wycombe (UK)
- RIKILT, Wageningen (NL)
- Royal Veterinary and Agricultural University, Copenhagen (DK)
- State Laboratory, Dublin (IRL)
- Technische Universität München, Bayerische Hauptversuchsanstalt für Landwirtschaft, Freising (D)
- TNO Biotechnology and Chemistry Institute, Zeist (NL)

ANALYTICAL METHODS USED

The methods used for certification involved instrumental determination by high performance liquid chromatography using a variety of separating and detection conditions and two-dimensional TLC. The methods also varied in their initial extraction and clean-up procedures. Details of the methods used are given in the certification report.

LEGAL NOTICE

This document was prepared under the sponsorship of the Commission of the European Communities. Neither the Commission of the European Communities nor any person acting on their behalf makes any warranty of representation, express or implied, that the use of any information, material, apparatus, method or process disclosed in this document may not infringe privately owned rights, or assume any liability with respect to the use of, or for damages resulting from the use of any information, material, apparatus, methods or process disclosed in this document.

NOTE

A detailed technical report on the analysis procedure and the treatment of the analytical data is supplied with each sample.

Certificate

Commission of the European Communities Community Bureau of Reference

bcr information

REFERENCE MATERIALS

THE CERTIFICATION OF THE AFLATOXIN B, MASS FRACTION OF THREE PEANUT MEAL AND TWO COMPOUND FEED REFERENCE MATERIALS

(CRMs 262, 263, 264, 375 & 376)

H.P. VAN EGMOND('), G. WOOD('), S. PATEL('), L. TUINSTRA('), A. BOENKE('), B. SCHURER('), P.J. WAGSTAFFE(')

Bilthoven
The Netherlands

(*) Leatherhead Food Research Association Leatherhead United Kingdom

(*) RHM Research and Engineering Ltd High Wycombe Bucks United Kingdom

(*) RIKILT Wageningen The Netherlands

(*) Commission of the European Communities Community Bureau of Reference 200 rue de la Loi B-1049 Brussels

Contract Nos 5394/1/5/227/90/9-BCR-NL (10) & 2721/1/5/227/87/5-BCR-NL (10)

FINAL REPORT

Directorate-General Science, Research and Development

EUR 15286 EN

Certification Report

Length: here 90 pages

Index of content: 4 pages

Content: feasibility study in terms of homogenity, stability,results, statistical evaluation, Certification ...

Why: the whole process has got to be traceable



1994

Matrix CRMs for Environmental Analysis

RTC is world leader in the production and certification of soil and sediment CRMs

Cooperation with US EPA

Available CRMs:

- Soil Matrix CRMs (Metals/Inorganic and Organic & Hydrocarbon)
- Soil, Sludge and Sediment Matrix CRMs Trace Element
- Underground Storage Tank Hydrocarbon Calibration CRMs
- Water Pollution Matrix (Metals/Inorganic and Organic & Hydrocarbon)
- And many more ...



Matrix CRMs for Environmental Analysis

Water CRMs

Matrix	Metals and Inorganics	Organics	Physical Properties
Water Supply (WS)	>35 products	>28 products	>18 products
Water Pollution (WP)	>29 products	>48 products	>20 products
Sea water	>3 products	>3 products	>5 products

Soil CRMs (Generic soils with known levels of analytes)

	Metals and Inorganics	Organics	Physical Properties
Soil / Solid waste	>13 products	>49 products	>7 products

Matrix CRMs (Soil with known sand, silt, clay and organic matter composition, and analytes)

	Metals and Inorganics	Organics	Physical Properties
Sandy Clay	✓	✓	✓
Loamy Sand	✓	✓	✓
Sandy Loam	✓	✓	✓
Sand	✓	✓	✓
Clay	✓	✓	✓
Sediment	✓	✓	_
Soil	✓	✓	_
Sandy Clay Loam	_	✓	_
Silty Clay	✓	✓	_
Loamy Clay	✓	✓	_
Clay Loam	✓	✓	_
Silty Loam	✓	✓	_
Loam	✓	✓	_
Clean Soil	✓	✓	✓
Other Matrices*	✓	✓	_

^{*}Other matrices include Fly Ash, (sewage) Sludge, Dusts and Paints. NB Paints are dry powders, and primarily certified for lead



Matrix CRMs for Environmental Analysis - CoA

Certificate of Am

Dioxins/Furans - Clay Loam

Number CRM981-10G Lot 019661 Bolvent (Matrix) Clay Loam

Hazard Inflant
Storage & Handling Store at room temperature.
Expiration Date: See Sample Label
Certification Date: December 04, 2012

SUPELCO

Solutions within.

Certified By: Christopher Rucinski - QA Director

Analyte	Unite	Certified 1,4 Value	RS	Stand. Deviat
1,2,3,4,6,7,8-Hpodf	P9/g	170 ± 6.82	2.05	15.7
1,2,3,4,7,8,9-Hpodf	pg/g	355 ± 10.7	2.05	23.8
1,2,3,4,6,7,8-Hpodd	P9/g	670 ± 31.2	2.05	72.7
Total Heptachlorodibenzo-p-dioxin (Total HPCDD)	pg/g	672 ± 32.0	2.05	72.9
Total Heptachlorodibenzofuran (Total HPCDF)	pg/g	526 ± 15.7	2.05	34.0
1,2,3,4,7,8-Hxcdd	Pg/g	51.7 ± 2.35	2.05	5.47
1,2,3,6,7,8-Hxedd	pg/g	509 ± 25.4	2.05	59.4
1,2,3,7,8,9-Hxcdd	pg/g	846 ± 46.3	2.05	109
Hoodd, total	P9/g	$1,390 \pm 69.7$	2.05	159
1,2,3,4,7,8-Hxcdf	P9/g	719 ± 32.0	2.05	74.2
1,2,3,6,7,8-Hxedf	pg/g	682 ± 32.7	2.05	76.3
1,2,3,7,8,9-Hxcdf	P9/g	302 ± 11.3	2.05	25.9
2,3,4,6,7,8-Headf	P9/g	564 ± 23.3	2.05	53.8
Total Hexachlorodibenzofuran (Total HxCDF)	pg/g	2,260 ± 92.0	2.05	207
1,2,3,4,6,7,8,9-OCDF	pg/g	190 ± 9.84	2.05	22.9
1,2,3,4,6,7,8,9-OCDD	P9/g	516 ± 29.3	2.05	69.3
1,2,3,7,8-Pecdd	P9/g	880 ± 30.3	2.05	69.3
1,2,3,7,8-Peodf	pg/g	794 ± 30.8	2.05	70.7
2,3,4,7,8-Peodf	P9/g	700 ± 24.5	2.05	55.9
Pecdf, total	P9/g	$1,510 \pm 47.9$	2.05	108
Peodd, total	P9/g	878 ± 31.6	2.05	71.1
TCDD, total	P9/9	108 ± 4.30	2.05	9.60

3

Analyte	Unite	Certified 1.4 Value	Æ	Steni Devis
2,3,7,8-TCDF	P9/g	651 ± 19.5	2.05	43.
TCDF, total	pg/g	660 ± 18.9	2.05	41.
2,3,7,8-TCDD	P9/g	107 ± 4.47	2.05	10:
PCDF, total	Pg/g	$5,240 \pm 226$	2.09	40
PCDD, total	pg/g	$3,640 \pm 214$	2.09	39

Additional Information

Description

The sample size provided is 10 g of soil.

The soil has been sterilized to minimize degradation of the sample.

The sample has been sized to 100 mesh.

Recommended storage condition is 4°C.

The sample has been intentionally prepared with an apparent headspace.

Storage

The sample should be stored at 4°C. It has been determined to be stable for the duration of After sub-sampling replace cap securely and store remaining sample at 4°C.

The shelf life of the product was determined by historic stability of similar CRM's. The explipional production is a stability to some stability of similar CRM's. The explipional production is a stability testing by a 17025 accredited laboratory.

Stability and shelf life after opening must be determined by the user, taking into account as

Preparation Instructions

Mix well prior to use.

For Method 8280A, an eight gram sample size contains sufficient concentrations of analyte quantitative limits will be realized.

For Method 8290, a one gram sample will contain concentrations of analytes such that levrealized. This holds true though levels of detection will be elevated owing to the reduced a

Scope and Application

The Dioxina and Funan in Soil Certified Reference Material (CRM) consists of a single poly of soil, furtified with dioxina and funans, being a natural matthe waste sample fie searsjet is interference, etc., so is typical for similar matrices necelved by the laboratory for analysis. I interference, set as the sample for similar matrices necelved by the laboratory for analysis, analysis that are listed on the enclosed Certificate of Analysis. The sample has been away round-orbit to meet the requirements specified by the ISO Guides 34 and 35, and ISO 176.

Evaluation of Results

The Reference Value, 95% confidence interval(CL), for the Reference Value and 95% Preobtained by the methics identified in the Scope and Application's section of this Certificate
from the beginning to the end of the bottling sequence and sent for analysis by an indepen
round-robin was used to calculate reference values by the USEPA EMBL-CINN's compute
The generated RNHEGHT mass. IMMEGHT standard deviation and IMMEGHT standard
Confidence Interval (CI) for the mean and the 95% Prediction Interval (PI). For normally dis
the classical calculation method used to generate a 95% CI. For non-Gaussian date sets I.
BIWEIGHT data are also used to calculate a 95% PI. The 95% PI compares well to a 95%
normally distributed data, the BIWEIGHT 95% PI typically represents approximately a 25 at
BIWEIGHT men. Again, the BIWEIGHT 95% PI typically represents expressionately a 25 at
Laboratories performing the same analytical procedures on a sample whose values have to
that the true mean, as determined by the method, is within the 95% CI window. Laboratorie
95% PI window 19 out of 20 analyses. Laboratories should use the PI as guidance for labo
Additional Information on the program may be obtained by referring to the reference or by
Additionally contact RTC for evaluations.





Page 2 of 3

Certificate of Analysis CERTIFIED REFERENCE MATERIAL

GERTIFIED REFERENCE MA

Dioxins/Furans - Clay Loam

Lot 019661 Solvent (Matrix) Clay Loam

Hazard Inftant

Storage & Handling Store at room temperature. Expiration Date: See Sample Label

Expiration Date See Sample Label Certification Date: December 04, 2012

Certified By: Christopher Rucinski - QA Director

Health and Safet Information

All RTC Certified Reference Materials are intended only for professional use by properly trained laboratory personnel. This CRM has been reviewed for both health and safety and shipping risks. It is classified as non hazardous and is not classified as hazardous goods for shipping by road, see or all hazardous.

A full international MSDS as a downloadable odf file is available at www.rt-corp.com

- 1 Certified values are the robust statistical mean when prepared according to instructions from an interfaboratory Study and internal rigorous teating. 2 The standard deviation is the robust statistical standard deviation from the round robin interfaboratory study.
- 2 The standard deviation is the robust statistical standard deviation from the round robin interiaconatory study.
 4 Expanded Uncertainty (Ucrm) All uncertainty values in this document expressed as ± value are expanded uncertainties.

5 A: Coverage factor derived from a t-distribution table, based on the degrees of freedom of the data set. Confidence interval = 95%

TRACEABILITY: The standard was manufactured under an ISO 17025 certified quality system. The balance used to weigh new materials is accurate to */0.0001g and calibrated regularly using mass standards traceable to NIST. All dilutions were preformed gravimetrically. Additionally, individual analytes are traceable to NIST SPMs where available and specified above.

HOMOGENEITY ASSESSMENT: Between-bottle homogeneity was assessed in accordance with 190 Guide 35. Completed units were sampled over the course of the bottling operation. Samples were taken in the following manner: the units produced in the bottling operation were divided into three chronological groups, those from the Early third, he Middle third, and the Late third (Groups). A pre-determined number of sample units were then randomly selected from each group. A subset of each group was then randomly selected for chemical snakysis. The results of the chemical snakysis were then compared by Single Factor Analysis of Verkinece (ANOVA).

UNCERTAINTY STATEMENT: Uncertainty values in this document are expressed as Expanded Uncertainty (Ucrm) corresponding to the 95% confidence interval. Ucrm is derived from the combined standard uncertainty multiplied by the coverage factor it, which is obtained from a f-distribution and degrees of feedom. The components of combined standard uncertainty include the uncertainties due to characterization, homogeneity, long term stability, and short term stability (transport). The components due to stability are generally considered to be negligible unless otherwise indicated by stability studies.

THIS PRODUCT WAS DESIGNED, PRODUCED AND VERIFIED FOR ACCURACY AND STABILITY IN ACCORDANCE WITH ISO 17025 (ACIssas Cert AT-1467) and ISO GUIDE 34 (ACIssas Cert AR-1470).

MSDG reports for components comprising greater than 1.0% of the solution or 0.1% for components known to be carcinogens are swallable upon request

Manufactured and certified by Sigma-Aldrich RTC, Inc.

318 - 3



Page 3 of 3

Matrix CRMs for Environmental Analysis - CoA

Certified Value with Uncertainty and Confidence Interval Additional Information

- Description
- Storage
- Preparation Instructions
- Scope and Application
- Evaluation of Results

Tracebility

Homogeneity Assessment

Uncertainty Assessment



Summary

New Method for Content Assignment of Organic Compounds by qNMR

- Traceability independent of chemical identity
- Low Uncertainties
- Non-destructive
- Structural information and information on impurities

Generation of new CRMs for Chromatography

- Quality Assurance
 - Ensuring Accuracy of the Mean and Precision
 - Ensuring an Analytical Method provides correct Results

Differentiation between pure substance CRMs/matrix CRMs



Thank you!



www.sigma-aldrich.com/organiccrm www.sigmaaldrich.com/rtccrm www.sigmaaldrich.com/ecrm



Lit. Code: QZN)