

# Tips and Tricks for Environmental GC Analysis

Dr. Frank Michel, Katherine K. Stenerson



[sigma-aldrich.com](http://sigma-aldrich.com)

# Troubleshooting Strategy


## • Approaching the problem...

- **Stop**, take a breath and think!
  - When did the problem start?
  - Has something changed?
- Check first to see if a “fix” for the problem is already known.
  - Talk to others in your lab
  - Check instrument maintenance/service log




# Troubleshooting Strategy

Isolate the source of the problem:



## First thing, review your method parameters...


- Injection type?
  - Should it be splitless?
- Split vent time
  - Too long or too short?
- Column flow
  - Using EPC?
- Heated zones
  - Double check temp. settings
- Liner type
  - Is there something better out there?



# Splitless Injection & Liners

# Splitless Injection & Liners


What you need to know about **splitless Injection**



## For splitless injections, consider:

- The splitless time  $\xrightarrow{\text{Too short}}$  Loss of response (esp. higher MW)  
 $\xrightarrow{\text{Too long}}$  Too much solvent on column

- The volume of the liner:



	B.P. (° C)	200° C Inlet Temp.		300° C Inlet Temp.	
		10 psi	30 psi	10 psi	30 psi
Methylene chloride	40	360	200	437	241
Methanol	64.5	570	315	691	382
Water	100	1279	706	1548	855

Resulting vapor volume of solvent

## Some Popular Styles for Splitless Injection

**Focus liner w/taper**



**2 mm ID, straight**



**Dual-tapered**



**Single-tapered**




**Semi-volatile  
Analysis**



# Liner ID

The ID of the liner can affect sensitivity:

The Use of a 2mm ID Liner will Increase Sensitivity for the Lighter Analytes




**Splitless injection,  
2mm vs. 4mm ID liner**

## Liner Care

If you *must* clean a liner....


- Handle with gloves or forceps.
- Use clean compressed gas and/or a fine brush to remove particles.
- Rinse in an appropriate solvent and dry with clean compressed gas.
- Use mineral acid and/or detergent only if absolutely necessary. Be sure to deactivate the liner after this process.
- If repacking with glass wool, make sure it has been deactivated.

### Undeactivated glass wool



1. 4,4'-DDE
2. Endrin
3. 4,4'-DDD
4. Endrin aldehyde
5. 4,4'-DDT
6. Endrin ketone


### Deactivated glass wool




10

## Give a new liner a chance!

Allow a “junk” injection after installing a new liner.



Freshly installed liner,  
first injection of  
pesticide standard  
(MSD)



Second injection of  
pesticide standard on  
same system

# Common Chromatographic Problems

# Common Chromatographic problems

## 1. Baseline Noise and Drift

Common causes:

- Column bleed
- Septa bleed
- Dirty detector
- Contaminants in carrier gas / carrier gas purity

## 2. Peak Shape & Response

- No response or poor response
- Extraneous peaks
- Poor peak shape




## Problem: Column Bleed

Did you know?

- Bleed results from the normal degradation of the stationary phase.
- All columns bleed to some extent.
- Bleed increases with temperature.
- The amount of bleed will increase in the presence of oxygen.

A Typical Bleed Profile:




## Column Bleed and an MSD


- Visible as baseline rise in the TIC.
- Check mass spectra for key bleed ions:
  - Stationary Phase -1: 73, 207, 281
  - Stationary Phase -5: 207, 281
  - Stationary Phase -1701: 207, 269
  - Stationary Phase -624: 207, 269

Make sure interface temp.  
< column max. temp.


### Common bleed ions



**D3**  
**D3 – CH<sub>3</sub>:**  
**207**



**D4**  
**D4 – CH<sub>3</sub>:**  
**281**



**TMS**  
**73**

## So, what can I do about bleed?

- Sufficiently purge column with carrier gas before ramping it up in temperature.
- Make sure carrier gas is scrubbed for water and oxygen.
- Check integrity of all fittings leading to the column.
- Do not heat the column above its maximum temp.
- Precondition the column prior to use.
- Use a high quality, high temperature septa and ferrules.



## To help prevent column bleed and other problems...remember gas purification

### Minimum recommendations for removal

	Carrier	Hydrogen	Air	Nitrogen	P-5
Oxygen	X				X
Water	X		X		
Hydrocarbons	X	X	X	X	X
Halocarbons					X



A wide variety of purifiers are available


## Problem: Too many peaks or “Ghost” Peaks




### Possible causes:

- Residue in the inlet liner and at the head of the column
- Contaminated syringe / and or wash solutions on an autosampler
- Sample carryover
- Contaminated carrier gas
- Septa pieces in liner

# Ghost peaks caused by septa pieces in a liner:




before



after

## Septa Bleed: MS Spectra



## Problem: Missing Peaks and Poor Response

Possible causes:

### Sample decomposition

- Activity in the inlet or column
- Injection port temperature too high
- Sample not stable enough for GC
- Standards not stable

### Column Installation

- Make sure your GC column is installed at the proper distance; both injector and detector


### Coelution

- Insufficient run time / final temperature
- Sample not volatile enough for GC
- Improper column installation




## Loss in response caused by creation of active sites:

Nasty samples can damage a column by creating active sites.




Responses of some acidic compounds were affected.





22

# Improper column installation can affect response and peak shape



Installation distance can affect response!



## What if my problem is coelution?

Go back to the basics – consider the resolution equation

$$R_s = \underbrace{(k/1+k)}_{\text{capacity}} \underbrace{(\alpha-1/\alpha)}_{\text{selectivity}} \underbrace{(N^{1/2}/4)}_{\text{efficiency}}$$

<b>column</b>	Longer column	→	capacity, efficiency
	Thicker film column	→	capacity
	Smaller ID column	→	efficiency
	Diff. Stationary phase	→	selectivity
<b>GC Parameters</b>	Carrier gas flow	→	efficiency
	Oven ramp rate	→	capacity
	Starting or ending oven temp.	→	capacity



## The best way to solve problems is to prevent them!

- Gas purification
  - Install and use appropriate filters/getters
- Injector maintenance
  - Liner, septa, seal
- Column installation
  - Check insertion distance
- Guard column
  - Use when necessary



# Troubleshooting Best Practices

## ✓ Documentation


- Use maintenance log books
- May save weeks of troubleshooting

## ✓ Make one change at a time

- To uncover root cause
- Multiple changes may offset each other

## ✓ Keep a 'good' trap

- Remove and store the trap as a reference for when issues occur at a later time
- Replace the caps, place in original shipping container, label properly, and protect from vibration



## Suggested Literature from Supelco

1. GC Column Selection Guide – Achieve Optimal Performance, T407133
2. Fast GC – A Practical Guide for Increasing Sample Throughput without Sacrificing Quality, T407096.
3. Capillary GC Inlet Liner Selection Guide (Bulletin 899A), T100899A
4. **Capillary GC Troubleshooting Guide: How to Locate Problems and Solve Them** (Bulletin 853C), T112853.
5. Purge and Trap System Guide (Bulletin 916), T197916
6. Gas Chromatography Accessories and Gas Purification/Management Products, T407103



## Capillary GC Troubleshooting Guide: How to Locate Problems and Solve Them

*The real task in correcting a problem with your capillary GC system is identifying the cause of the problem without wasting time. The systematic approach to troubleshooting described in this guide will enable you to solve many problems yourself. The guide also contains suggestions for maintaining your system, including the column, at optimal performance levels. By following these recommendations, you can reduce repair costs and instrument down time.*



### Contents

☐ Troubleshooting Suggestions	2
☐ Isolating the Problem Source	2-6
☐ Tips to Help Prevent Problems	6-8
☐ Troubleshooting Table	9-28
☐ Appendix	29, 30
☐ Helpful Products	31-40

Symptom	Possible Cause	Remedy
<b>Symptom No. 2 — Missing peaks/solvent peak only (cont'd.)</b>		
 <p>796-0217</p>	5. Incorrect linear velocity/column flow rate. 6. Sample components adsorbed by column or inlet liner.  7. Column cannot separate components from solvent.	5. Measure column flow and adjust if necessary (see page 4). 6. Inject standard on column known to be performing well. If results are good, rinse original column (bonded phase only) or remove 1-2 coils from inlet end (see page 8). If column performance is not restored, replace column. To prolong column life, use guard columns (see page 8) and reverse column periodically. Remove inlet liner and check cleanliness. Replace glass wool and packing or use new, deactivated liner. If sample has never been analyzed and is chemically active, you may need a column with a different (or specialized) stationary phase. 7. Change solvent or column.
<b>Symptom No. 3 — Detector response low for all peaks (retention times correct)</b>		
 <p>796-0218</p>	1. Sample size too small or split ratio too large. 2. Sensitivity setting wrong. 3. Makeup gas flow inadequate. 4. Poor injection technique. 5. Syringe defective. 6. Carrier gas leak at septum or column connection. 7. Injection port temperature too low (sample not vaporized). 8. FID: hydrogen or air flow incorrect. 9. FID: low oxygen level in compressed air.	1. Increase sample size or reduce split ratio. 2. Check sensitivity setting, adjust if necessary. Inject standard for comparison. 3. Adjust flow to detector manufacturer's specifications. 4. Use correct syringe size and appropriate injection technique for sample and analysis. 5. Inject sample with a new syringe. If problem disappears, discard old syringe. 6. Check for leaks (see page 3). Replace septum or tighten connections if necessary. 7. Verify temperature with accurate thermometer, adjust if necessary. If instrument setting and thermometer agree, increase temperature (do not exceed stationary phase limit). 8. Adjust gas flows to detector manufacturer's specifications. 9. Replace air tank.

**Thank you!**

