

# Application of Solid Phase Micro Extraction in Environmental Analysis

Frank Michel



[sigma-aldrich.com/analytical](http://sigma-aldrich.com/analytical)



## SPME in Environmental Analysis

More than 900 references for Environmental Analysis on currently edition of SPME CD:

- Water
- Soil
- Air
- Pesticides

Plus literature, applications and videos





## Official Methods in Environmental Analysis using SPME

ASTM D 6520, 2000

- Standard Practice for the SPME of Water and its Headspace for the Analysis of Volatile and Semi-Volatile Organic Compounds

ASTM D 6889, 2003

- Standard Practice for Fast Screening for Volatile Organic Compounds in Water Using Solid Phase Microextraction (SPME)

EPA Method 8272, 2007

- Parent and Alkyl PAHs in Sediment Pore Water by SPME-GC/MS

ISO 27108 (DIN 38407-34), 2013

- Determination of selected plant treatment agents and biocide products - Method using SPME followed by GC-MS

ISO 17943 (DIN 38407-41)

- Determination of VOCs in water - GC-MS after HS-SPME



## ISO 27108 (DIN 38407-34)

Determination of selected plant treatment agents and biocide products in drinking water, ground water and surface water - Method using SPME followed by GC-MS

- Dichlobenil
- Desethylatrazin
- Desethylterbutylazin
- Simazin
- Atrazin
- Lindan
- Terbutylazin
- Metribuzin
- Parathion-methyl
- Heptachlor
- Terbutryn
- Aldrin
- Metolachlor
- Parathion-ethyl
- exo-Heptachlorepoxyd
- Pendimethalin
- endo-Heptachlorepoxyd
- Triclosan
- Dieldrin
- Carfentrazone-ethyl
- Diflufenican
- Mefenpyr-diethyl
- Method might be suitable for other compounds – need to be evaluated individually
- Operational range of method: above 0.05 µg/L (depending on matrix)
- IS: Atrazin-d5 or Lindan-d6



## Conditions

### SPME Conditions

- Fiber: Polyacrylate
- Conditioning: Fiber 20 min @ 300° C
- Extraction: pH 6-8, NaCl close to saturation (e.g. 2.4 g/8 mL sample)  
60 min @ 30° C, stirring 250 rpm
- Desorption: 10 min @ 280° C , splitless (6 min)

### GC conditions (example)

- Injection: KAS/CIS 60° C; 10° C/s to 280° C; 10min; Splitless
- Column: 5% Phenyl, 30m x 0.25mm, 0.25µm
- Carrier gas: Helium (5.0), 0,9mL/min
- Oven: 60° C, 2 min; 4° C/min to 220° C; 10° C/min to 300° C
- Detector: MS, EI 70eV, SIM




## Selected Ions for Identification and Quantification

| Analytes              | Selected Ions for Identification and Quantification (m/z) |
|-----------------------|---|
| Dichlobenil           | 100, 136, 171, 173  |
| Desethylatrazin       | 145, 172, 174, 187  |
| Desethylterbutylazin  | 186, 188, 201   |
| Simazin               | 173, 186, 201   |
| Atrazin               | 173, 200, 215   |
| Lindan                | 109, 181, 183, 219  |
| Terbutylazin          | 173, 214, 229   |
| Metribuzin            | 103, 144, 198, 214  |
| Parathion-methyl      | 109, 125, 263   |
| Heptachlor            | 237, 272, 274, 337  |
| Terbutryn             | 170, 185, 226, 241  |
| Aldrin                | 261, 263, 265, 293  |
| Metolachlor           | 162, 238, 240   |
| Parathion-ethyl       | 109, 155, 263, 291  |
| exo-Heptachlorepoxyd  | 353, 355, 357   |
| Pendimethalin         | 162, 191, 252, 281  |
| endo-Heptachlorepoxyd | 183, 253, 289   |
| Triclosan             | 218, 288, 290   |
| Dieldrin              | 79, 263, 277, 279   |
| Carfentrazon-ethyl    | 290, 312, 340, 411  |
| Diflufenican          | 246, 266, 394   |
| Mefenpyr-diethyl      | 227, 253, 255, 299  |




## Example Chromatogram





## Example Chromatogram








## Peak Identification in Example Chromatogram


- |                        |                          |
|------------------------|--------------------------|
| 1 Dichlobenil          | 13 Aldrin                |
| 2 Desethylatrazin      | 14 Metolachlor           |
| 3 Desethylterbutylazin | 15 Parathion-ethyl       |
| 4 Simazin              | 16 exo-Heptachlorepoxid  |
| 5 Atrazin-d5           | 17 Pendimethalin         |
| 6 Atrazin              | 18 endo-Heptachlorepoxid |
| 7 Lindan               | 19 Triclosan             |
| 8 Terbutylazin         | 20 Dieldrin              |
| 9 Metribuzin           | 21 Carfentrazon-ethyl    |
| 10 Parathion-methyl    | 22 Diflufenican          |
| 11 Heptachlor          | 23 Mefenpyr-diethyl      |
| 12 Terbutryn           |                          |




## Challenge: highly polar Pesticides



Mecoprop (MCP)



Bromoxynil



Bentazon



## Challenge: highly polar Pesticides

Solution: Combination of SPME and Derivatization


- Derivatization of analytes in water sample *before the extraction*
- Derivatization of analytes after extraction
  - On fiber
  - In injector

From presentation: Dr. Friedrich Werres,  
„Pesticides in water – efficient determination by SPME-GC/MS“ (2005)




## Methylation for highly polar Pesticides

### Derivatization with Diazomethane




### Derivatization with TMSH



From presentation: Dr. Friedrich Werres,  
„Pesticides in water – efficient determination by SPME-GC/MS“ (2005)




## Methylation for highly polar Pesticides - TMSH



From presentation: Dr. Friedrich Werres,  
„Pesticides in water – efficient determination by SPME-GC/MS“ (2005)



## Methylation for highly polar Pesticides - Diazomethane



From presentation: Dr. Friedrich Werres,  
„Pesticides in water – efficient determination by SPME-GC/MS“ (2005)



## Nitrosamines in water

SPME used in a direct immersion extraction of nitrosamines from water. These are often difficult analytes to extract especially dimethylnitrosoamine.

- Nitrosodimethylamine
- Nitrosodiethylamine
- Nitrosomethylethylamine
- Nitrosodipropylamine
- Nitrosopiperidine
- Nitrosodibutylamine
- Nitrosodiphenylamine





## 10ppb Nitrosamines in Water: SPME-GC/MS


|             |  |
|-------------|--|
| Sample:     | analytes in (water + 25% KCl, pH 10)                 |
| SPME Fiber: | 65 $\mu$ m PDMS-DVB                                  |
| Extraction: | immersion, 15 min (rapid stirring)                   |
| Desorption: | 270° C, 1 min  |
| Column:     | PTA-5 (amine deactivated, 30m x 0.32mm, 0.5 $\mu$ m) |
| Oven:       | 50° C (1 min) to 250° C at 10° C/min, hold 2 min     |
| Carrier:    | helium, 30cm/sec                                     |
| Det.:       | GC/MS (quadrupole, SIM)                              |
| Inj.:       | splitless, 250° C (0.75mm ID liner)                  |





## 10ppb Nitrosamines in Water: SPME-GC/MS

1. Nitrosodimethylamine
2. Nitrosodiethylamine
3. Nitrosomethylethylamine
4. Nitrosodipropylamine
5. Nitrosopiperidine
6. Nitrosodibutylamine
7. Nitrosodiphenylamine



Chromatogram courtesy of J. Clark, Liggett Group, Inc.



## Standard Method 6040D: Odor compounds in water

This application demonstrates the use of headspace SPME and the SLB-5ms column for the low level extraction and detection of odor compounds in drinking water. Specifically, SPME is described for this analysis in Standard Method 6040D of these compounds at the part-per-trillion (ppt) level.





## Odor-causing compounds in water

### SPME

- SPME fiber: 2 cm Metal 50/30  $\mu\text{m}$  DVB/CAR/PDMS
- Extraction: headspace, 65 ° C (30 min.)
- Desorption: 3 min. at 260 ° C


### GC-MS

- Column: SLB-5ms, 30 m x 0.25 mm I.D., 0.25  $\mu\text{m}$
- Oven: 60 ° C (2 min.), 8 ° C/min. to 200 ° C
- MSD interf.: 300 ° C
- Scan range: SIM, m/z = 137, 124, 95, 112
- Carrier gas: helium, 1 mL/min. constant
- Liner: 0.75 mm I.D., SPME
- Sample: 20/10 ppt odor comps. in 25 mL water + 25 % NaCl



## Odor-causing compounds in water


1. 2-isopropyl-3-methoxypyrazine, 20 ppt
2. 2-isobutyl-3-methoxypyrazine, 20 ppt
3. 2-methylisoborneol, 10 ppt
4. (+/-)Geosmin, 10 ppt





## Odor-causing compounds in water at 2 ppt (GC/MS)

1. 2-Isopropyl-3-methoxypyrazine (IPMP)
2. 2-Isobutyl-3-methoxypyrazine (IBMP)
3. 2-Methylisoborneol (MIB)
4. 2,4,6-Trichloroanisole (I.S. 8ppt)
5. ( $\pm$ ) Geosmin





## ISO Standard 17943

Water quality - Determination of VOCs in water - Method using HS-SPME followed by GC-MS

- Broad use of VOCs (Volatile Organic Compounds) in many products, more than 60 are covered in this method
  - Halogenated hydrocarbons
  - Trihalogen methanes
  - Gasoline additives (like BTEX, MTBE and ETBE)
  - Naphthalene
  - 2-Ethyl-4-methyl-1,3-dioxolane and highly odorous substances like geosmin and 2-methylisoborneol
- Concern for human health as many of them are toxic and known or suspected to be carcinogenic
- Application in drinking water, ground water and surface water



## VOCs in Water

Water Framework Directive requires application of ISO or CEN standards (standard methods)

Current methods have not been state-of-the-art:

- ISO 10301 (1997, Liquid-Liquid extraction, GC/FID or GC/ECD)
- ISO 11423 (1997, Headspace-GC/FID or GC/ECD)
  - ☹ Sensitivity & Selectivity
- ISO 15680 (2003, Purge and Trap, GC/MS)
  - ☹ Susceptibility to contamination, automation is complex



## ISO 17943 – Content

- 45 pages
- Practical tips on using SPME and GC/MS
- Scope and principle
- Reagents, apparatus, sampling and sample pretreatment, procedure, calibration
- Calculation and expression of the results
- Examples of suitable SPME fibres, GC columns, internal standards, gas chromatographic conditions and example chromatograms






## SPME Method (Supelco Applications Lab)

|                             |  |
|-----------------------------|--|
| <u>Sample volume:</u>       | 10 mL  |
| <u>HS-Vial:</u>             | 20 mL, addition of 3 g salt  |
| <u>SPME fiber:</u>          | DVB/CAR/PDMS, 24 gauge   |
| <u>Incubation time:</u>     | 10 min @ 40 ° C  |
| <u>Extraction time:</u>     | 10 min @ 40 ° C  |
| <u>Autosampler:</u>         | CTC Combi PAL (agitated by circular motion of the vial, velocity: 250 rpm) |
| <u>Desorption/Injector:</u> | 10 min @ 270 ° C   |



# GC Method (Supelco Applications Lab)



**GC:** Varian CP-3800  
**Column:** VOCOL, 60m x 1,5 $\mu$ m x 0,25mm  
**Carrier gas:** He  
**Flow:** 1 mL/min  
**Injection:** Splitless, SPME liner w/ 0,75 mm ID  
**Oven program:** 35 C, 1 min; 10 C/min to 150 C; 20 C/min to 250, 20 min  
**Sample:** 61 VOCs, 1 ppm, in water plus three internal standards



## Interlaboratory Trial for Validation: Participants

- Austria 1x
- Brazil 2x
- Canada 2x
- Croatia 1x
- France 2x
- Germany 12x
- Great Britain 1x
- Italy 5x
- Portugal 3x
- Romania 1x
- Serbia 1x
- South Africa 2x
- Spain 4x
- Sweden 1x
- Swizerland 1x
- United States 3x

**42 Participants out of 16 countries**





## Interlaboratory Trial

Determination of the concentration of 61 compounds in the two samples

Four independent replicate analysis from each of the 2 samples

Strictly following the procedure as prescribed in the draft standard (ISO/CD 17943)

Results had to be delivered 30 days after receipt of the samples





## Samples for Interlaboratory Trial

**Sample 1: Surface water** was taken from an urban and industrialized area (river Ruhr in Muelheim, Germany)

- Filtration using a glass fiber filter.
- Stabilization with 50 mg/L sodium azide.

**Sample 2: Municipal wastewater** was taken from a plant effluent.

- Sedimentation and pumping into a large fluid tank while being filtrated by both 5  $\mu\text{m}$  and 1  $\mu\text{m}$  and irradiated by UV
- Sterilisation (80 ° C), introduction of gas: (1) CO<sub>2</sub>, (2) N<sub>2</sub>
- Stabilization with 50 mg/L sodium azide.

### Spiking of samples:

- Surface water: 0,02 – 0,80  $\mu\text{g/l}$  (~ 50 % < 0,10 g/l)
- Waste water: 0,05 – 3,00  $\mu\text{g/l}$  (~ 50 % < 0,50 g/l)

Samples were tested for homogeneity and stability





## Results from Interlaboratory Trial

No submission of results: 9 labs

Significant deviation from the procedure prescribed: 6 labs

- calibration without internal standards (3x)
- other major deviations from draft ISO/CD 17943 (3x)

**A total of 27 labs reported results to be included in the evaluation process according ISO 5725-2**


- All parameters analysed: 10 labs
- Nearly all parameters analysed: 9 labs
- Nearly each parameter had been analysed by > 20 labs



# Results from Interlaboratory Trial

Probe: 1  
Merkmal: 2-chlorotoluene  
Referenzwert: 0,100 µg/l

Mittelwert: 0,096 µg/l  
Vergleich-Stdabw. (SR): 0,022 µg/l  
Wiederhol-Stdabw. (Sr): 0,005 µg/l



PROLab Smart für ISO 5725-2



## Results from Interlaboratory Trial

Analysis for

- Recovery rate (from assigned value)
  - For most of the compounds between 84 and 116 % (surface water) and 81 and 118 % (waste water)
- Reproducibility standard deviation
  - For most of the compounds less than 31 % (surface water) and less than 35 % (waste water)
- Repeatability standard deviation
  - For most of the compounds less than 10 % (surface water) and less than 8 % (waste water)





Environmental

**Thank you!**

**SPME**  
  
**SUPELCO**

**S ample**  
**P rep**  
**M ade**  
**E asy**